



BPAP 30



BPAP A30

(Automatic)  
Bi-level Non-invasive Ventilator

Sepray<sup>®</sup> BPAP 30  
BPAP A30

## Healthy & Self-confidence


Breathing & Sleeping  
Solution for OSAHS patients




RESPITEK

Quality • Service


Therapy Interface



Info. Interface



Real-time display treatment data  
3.5" color display



## Design

Two-in-one structure 

Main unit and humidifier are combined into one structure

Compact and Portable

Dimensions: 24.3(L) x 16.7(W) x 13(H) cm

Weight: 1.65kg



## Special Humidifier

Patented 360° anti-countercurrent water tank design

Humidifier temperatures adjustable





## Function

COMF, Pressure Release Technology

Ensure treatment comfort degree, improve treatment compliance.

Blood oxygen monitoring function

Simultaneous display of blood oxygen monitoring data.

Multi-Parameter display

Real-time monitoring of main professional treatment data.

High-performance blower

Stable, quiet and efficient air leak compensation.

Scientific treatment algorithm

Ability to monitor and identify a variety of respiratory events automatically and accurately, response to apnea, hypopnea, snore, large leaks, continuous airflow limitation etc. events immediately.

| Parameter/Model  | BPAP 30                                  | BPAP A30                                     |
|------------------|--|--|
| Mode             | CPAP<br>Bilevel                          | CPAP<br>Auto CPAP<br>Bilevel<br>Auto Bilevel |
| Pressure         | 4cmH <sub>2</sub> O–30cmH <sub>2</sub> O | 4cmH <sub>2</sub> O–30cmH <sub>2</sub> O     |
| Max PS           | no limits                                | Auto Bi-level less than 8cmH <sub>2</sub> O  |
| COMF             | 1–3 levels                               | 1–3 levels                                   |
| Rise time        | 1–4 levels ( 200ms–500ms )               | 1–4 levels ( 200ms–500ms )                   |
| Ramp time        | 0–45min adjustable                       | 0–45min adjustable                           |
| Humidifier       | 1–5 levels                               | 1–5 levels                                   |
| Auto ON/OFF      | Available                                | Available                                    |
| Alert            | Available                                | Available                                    |
| LED              | Available                                | Available                                    |
| WiFi             | Available                                | Available                                    |
| Data storage     | TF card (4G)                             | TF card (4G)                                 |
| Split Night      | Unavailable                              | Available ( 0–240min )                       |
| Language         | CN/EN                                    | CN/EN  |
| SPO <sub>2</sub> | Optional configuration                   | Optional configuration                       |



**İstanbul Merkez Ofis**  
 Altunizade Mah. Kısıklı Cad. No:35/1 Kat:1 Üsküdar / İstanbul  
 Tel : 0216 545 8080

[www.respitek.com.tr](http://www.respitek.com.tr)



**Micome Medical Technology Development Co., Ltd.**

E-mail: [info@micome.cc](mailto:info@micome.cc) Tel:+86-0731-89912333-877 Fax: +86-0731-89916333 Website: [www.micome.cc](http://www.micome.cc)  
 Add: No.8,Lujing Road,High-Tech Zone,Yuelu,Changsha,China